

TTBK2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200 IF 1:20 - 1:100

Calculated MW: 137kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human TTBK2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TTBK; SCA11;

Catalog #: A5168

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 146057

Isotype: IgG

Swiss Prot: Q6IQ55

Purity: Affinity purification

For research use only.

Background:

This gene encodes a serine-threonine kinase that putatively phosphorylates tau and tubulin proteins. Mutations in this gene cause spinocerebellar ataxia type 11 (SCA11); a neurodegenerative disease characterized by progressive ataxia and atrophy of the cerebellum and brainstem.

To place an order, please [Click HERE](#).