

COMMD1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

Calculated MW: 21kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human COMMD1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MURR1; C2orf5;

Catalog #: A5174

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 150684

Isotype: IgG

Swiss Prot: Q8N668

Purity: Affinity purification

For research use only.

Background:

COMMD1 is a regulator of copper homeostasis, sodium uptake, and NF-kappa-B (see MIM 164011) signaling (de Bie et al., 2005)

To place an order, please [Click HERE](#).