

FMR1NB

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:29kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human FMR1NB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CT37; NYSAR35; NY-SAR-35;

To place an order, please [Click HERE](#).

Catalog #:A5183

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:158521

Isotype:IgG

Swiss Prot:Q8N0W7

Purity:Affinity purification

For research use only.