

SLFN1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 46kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human SLFN1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A5198

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 200172

Isotype: IgG

Swiss Prot: Q499Z3

Purity: Affinity purification

For research use only.