

ARL13B

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 49kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ARL13B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

JBTS8; ARL2L1;

Catalog #: A5200

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 200894

Isotype: IgG

Swiss Prot: Q3SXY8

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the ADP-ribosylation factor-like family. The encoded protein is a small GTPase that contains both N-terminal and C-terminal guanine nucleotide-binding motifs.

This protein is localized in the cilia and plays a role in cilia formation and in maintenance of cilia.

Mutations in this gene are the cause of Joubert syndrome 8. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).