

USP12

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 43kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human USP12

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

UBH1; USP12L1;

To place an order, please [Click HERE](#).

Catalog #: A5201

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 219333

Isotype: IgG

Swiss Prot: O75317

Purity: Affinity purification

For research use only.