

TMED4

Reactivity: Human Mouse Rat

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:10 - 1:100

Calculated MW: 26kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human TMED4

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HNLF; ERS25; GMP25iso; p24alpha3;

To place an order, please [Click HERE](#).

Catalog #: A5208

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 222068

Isotype: IgG

Swiss Prot: Q7Z7H5

Purity: Affinity purification

For research use only.