

LCORL

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 67kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human LCORL

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MLR1;

Catalog #: A5209

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 254251

Isotype: IgG

Swiss Prot: Q8N3X6

Purity: Affinity purification

For research use only.

Background:

This gene encodes a transcription factor that appears to function in spermatogenesis.

Polymorphisms in this gene are associated with measures of skeletal frame size and adult height.

Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).