

## MMAB

**Reactivity:** Human Mouse

**Tested applications:** WB IHC IP

**Recommended Dilution:** WB 1:200 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

**Calculated MW:** 27kDa

**Observed MW:** Refer to figures

**Immunogen:**

A synthetic Peptide of human MMAB

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ATR; cob; cblB; CFAP23;

**Catalog #:** A5221

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 326625

**Isotype:** IgG

**Swiss Prot:** Q96EY8

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a protein that catalyzes the final step in the conversion of vitamin B(12) into adenosylcobalamin (AdoCbl), a vitamin B12-containing coenzyme for methylmalonyl-CoA mutase. Mutations in the gene are the cause of vitamin B12-dependent methylmalonic aciduria linked to the cblB complementation group. Alternatively spliced transcript variants have been found.

*To place an order, please [Click HERE](#).*