

TSPAN33

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:1000 IHC 1:50 - 1:200

Calculated MW: 32kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TSPAN33

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PEN;

To place an order, please [Click HERE](#).

Catalog #: A5222

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 340348

Isotype: IgG

Swiss Prot: Q86UF1

Purity: Affinity purification

For research use only.