

SMTNL2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 50kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human SMTNL2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A5224

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 342527

Isotype: IgG

Swiss Prot: Q2TAL5

Purity: Affinity purification

For research use only.