

GSTK1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 25kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GSTK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GST; GST13; hGSTK1; GSTK1-1; GST13-13; GST 13-13;

Catalog #: A5226

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 373156

Isotype: IgG

Swiss Prot: Q9Y2Q3

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the kappa class of the glutathione transferase superfamily of enzymes that function in cellular detoxification. The encoded protein is localized to the peroxisome and catalyzes the conjugation of glutathione to a wide range of hydrophobic substrates facilitating the removal of these compounds from cells. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).