

CISD2

Reactivity: Human Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 15kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CISD2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

ERIS; WFS2; ZCD2; NAF-1; Miner1;

Catalog #: A5231

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 493856

Isotype: IgG

Swiss Prot: Q8N5K1

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a zinc finger protein that localizes to the endoplasmic reticulum. The encoded protein binds an iron/sulfur cluster and may be involved in calcium homeostasis. Defects in this gene are a cause of Wolfram syndrome 2.

To place an order, please [Click HERE](#).