

CPLX3

Reactivity: Human Mouse Rat

Tested applications: WB IF IP

Recommended Dilution: WB 1:200 - 1:2000 IF 1:10 - 1:100 IP 1:20 - 1:50

Calculated MW: 18kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human CPLX3

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CPXIII; CPX-III; Nbla11589;

To place an order, please [Click HERE](#).

Catalog #: A5234

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 594855

Isotype: IgG

Swiss Prot: Q8WVH0

Purity: Affinity purification

For research use only.