

LGALS7B

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 15kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human LGALS7B

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PI7; GAL7; Gal-7; HKL-14; LGALS7;

Catalog #: A5235

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 653499

Isotype: IgG

Swiss Prot: P47929

Purity: Affinity purification

For research use only.

Background:

The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. Differential and in situ hybridization studies indicate that this lectin is specifically expressed in keratinocytes and found mainly in stratified squamous epithelium. A duplicate copy of this gene (GeneID:3963) is found adjacent to, but on the opposite strand on chromosome 19.

To place an order, please [Click HERE](#).