

TOMM6

Reactivity:Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:8kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TOMM6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

h

Synonym:

OBTP; TOM6;

To place an order, please [Click HERE](#).

Catalog #:A5240

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:100188893

Isotype:IgG

Swiss Prot:Q96B49

Purity:Affinity purification

For research use only.