

SNAI1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 29kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SNAI1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

SNA; SNAH; SNAIL; SLUGH2; SNAIL1; dJ710H13.1;

Catalog #: A5243

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6615

Isotype: IgG

Swiss Prot: O95863

Purity: Affinity purification

For research use only.

Background:

The Drosophila embryonic protein SNAI1 is a zinc finger transcriptional repressor which downregulates the expression of ectodermal genes within the mesoderm. The nuclear protein encoded by this gene is structurally similar to the Drosophila snail protein, and is also thought to be critical for mesoderm formation in the developing embryo. At least two variants of a similar processed pseudogene have been found on chromosome 2.

To place an order, please [Click HERE](#).