

ERAP1

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 107kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ERAP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

a

Synonym:

A-LAP; ALAP; APPILS; ARTS-1; ARTS1; ERAAP; ERAAP1; KIAA0525; PILS-AP; PILSAP;

Catalog #: A5248

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51752

Isotype: IgG

Swiss Prot: Q9NZ08

Purity: Affinity purification

For research use only.

Background:

ERAP1 is an aminopeptidase involved in trimming HLA class I-binding precursors so that they can be presented on MHC class I molecules. It acts as a monomer or as a heterodimer with ERAP2.

ERAP1 may also be involved in blood pressure regulation by inactivating angiotensin II.

To place an order, please [Click HERE](#).