

CD99L2

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW: 28kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human CD99L2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD99B; MIC2L1

Catalog #: A5288

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 83692

Isotype: IgG

Swiss Prot: Q8TCZ2

Purity: Affinity purification

For research use only.

Background:

This gene encodes a cell-surface protein that is similar to CD99. A similar protein in mouse functions as an adhesion molecule during leukocyte extravasation. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).