

CHRNA1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 31kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CHRNA1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

g

Synonym:

ACHRA; ACHRD; CHRNA; CMS2A; FCCMS; SCCMS

Catalog #: A5306

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1134

Isotype: IgG

Swiss Prot: P02708

Purity: Affinity purification

For research use only.

Background:

The muscle acetylcholine receptor consists of 5 subunits of 4 different types: 2 alpha subunits and 1 each of the beta, gamma, and delta subunits. This gene encodes an alpha subunit that plays a role in acetylcholine binding/channel gating. Alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).