

FABP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 14kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human FABP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FABPL; L-FABP

Catalog #: A5311

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2168

Isotype: IgG

Swiss Prot: P07148

Purity: Affinity purification

For research use only.

Background:

This gene encodes the fatty acid binding protein found in liver. Fatty acid binding proteins are a family of small, highly conserved, cytoplasmic proteins that bind long-chain fatty acids and other hydrophobic ligands. This protein and FABP6 (the ileal fatty acid binding protein) are also able to bind bile acids. It is thought that FABPs roles include fatty acid uptake, transport, and metabolism.

To place an order, please [Click HERE](#).