

GTF2A1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 19kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human GTF2A1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TF2A1; TFIIA; TFIIAL; TFIIA-42;

Catalog #: A5345

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2957

Isotype: IgG

Swiss Prot: P52655

Purity: Affinity purification

For research use only.

Background:

Accurate transcription initiation on TATA-containing class II genes involves the ordered assembly of RNA polymerase II (POLR2A; MIM 180660) and several general initiation factors (summarized by DeJong and Roeder, 1993 [PubMed 8224848]). One of these factors is TFIIA, which when purified from HeLa extracts consists of 35-, 19-, and 12-kD subunits.

To place an order, please [Click HERE](#).