

NEK2

Reactivity:Human Mouse

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:51kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human NEK2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

NLK1; NEK2A; HsPK21;

Catalog #:A5355

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4751

Isotype:IgG

Swiss Prot:P51955

Purity:Affinity purification

For research use only.

Background:

This gene encodes a serine/threonine-protein kinase that is involved in mitotic regulation. This protein is localized to the centrosome, and undetectable during G1 phase, but accumulates progressively throughout the S phase, reaching maximal levels in late G2 phase. Alternatively spliced transcript variants encoding different isoforms with distinct C-termini have been noted for this gene.

To place an order, please [Click HERE](#).