

ID3

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 13kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ID3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HEIR-1; bHLHb25;

Catalog #: A5375

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3399

Isotype: IgG

Swiss Prot: Q02535

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a helix-loop-helix (HLH) protein that can form heterodimers with other HLH proteins. However, the encoded protein lacks a basic DNA-binding domain and therefore inhibits the DNA binding of any HLH protein with which it interacts.

To place an order, please [Click HERE](#).