

AMPH

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 76kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human AMPH

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AMPH1

Catalog #: A5389

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 273

Isotype: IgG

Swiss Prot: P49418

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein associated with the cytoplasmic surface of synaptic vesicles. A subset of patients with stiff-man syndrome who were also affected by breast cancer are positive for autoantibodies against this protein. Alternate splicing of this gene results in two transcript variants encoding different isoforms. Additional splice variants have been described, but their full length sequences have not been determined. A pseudogene of this gene is found on chromosome 11.

To place an order, please [Click HERE](#).