

AP2A2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 104kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human AP2A2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

i

Synonym:

HIP9; HYPJ; ADTAB; HIP-9; CLAPA2

Catalog #: A5391

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 161

Isotype: IgG

Swiss Prot: O94973

Purity: Affinity purification

For research use only.

Background:

Adaptins are components of the adaptor complexes which link clathrin to receptors in coated vesicles. Clathrin-associated protein complexes are believed to interact with the cytoplasmic tails of membrane proteins, leading to their selection and concentration. Alpha adaptin is a subunit of the plasma membrane adaptor. It binds polyphosphoinositides.

To place an order, please [Click HERE](#).