

ASH2L

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:50

Calculated MW: 69kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ASH2L

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

ASH2; Bre2; ASH2L1; ASH2L2;

To place an order, please [Click HERE](#).

Catalog #: A5392

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9070

Isotype: IgG

Swiss Prot: Q9UBL3

Purity: Affinity purification

For research use only.