

PRKRA

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 34kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PRKRA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RAX; PACT; DYT16; HSD14;

Catalog #: A5417

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8575

Isotype: IgG

Swiss Prot: O75569

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein kinase activated by double-stranded RNA which mediates the effects of interferon in response to viral infection. Mutations in this gene have been associated with dystonia. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).