

SERPINB3

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 44kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SERPINB3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SCC; T4-A; SCCA1; SCCA-1; HsT1196; SCCA-PD;

To place an order, please [Click HERE](#).

Catalog #: A5418

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6317

Isotype: IgG

Swiss Prot: P29508

Purity: Affinity purification

For research use only.