

TFDP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 45kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TFDP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bk

Synonym:

DP1; Dp-1; DRTF1;

Catalog #: A5422

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7027

Isotype: IgG

Swiss Prot: Q14186

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of a family of transcription factors that heterodimerize with E2F proteins to enhance their DNA-binding activity and promote transcription from E2F target genes. The encoded protein functions as part of this complex to control the transcriptional activity of numerous genes involved in cell cycle progression from G1 to S phase. Alternative splicing results in multiple transcript variants. Pseudogenes of this gene are found on chromosomes 1, 15, and X.

To place an order, please [Click HERE](#).