

TRIB3

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 40kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TRIB3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

NIPK; SINK; TRB3; SKIP3; C20orf97;

Catalog #: A5424

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57761

Isotype: IgG

Swiss Prot: Q96RU7

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a putative protein kinase that is induced by the transcription factor NF-kappaB. The encoded protein is a negative regulator of NF-kappaB and can also sensitize cells to TNF- and TRAIL-induced apoptosis. In addition, this protein can negatively regulate the cell survival serine-threonine kinase AKT1. Differential promoter usage and alternate splicing result in multiple transcript variants.

To place an order, please [Click HERE](#).