

TMSB4X

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 5kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TMSB4X

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FX; TB4X; PTMB4; TMSB4;

Catalog #: A5438

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7114

Isotype: IgG

Swiss Prot: P62328

Purity: Affinity purification

For research use only.

Background:

This gene encodes an actin sequestering protein which plays a role in regulation of actin polymerization. The protein is also involved in cell proliferation, migration, and differentiation. This gene escapes X inactivation and has a homolog on chromosome Y.

To place an order, please [Click HERE](#).