

PI3

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 12kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PI3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Cementoin; Elafin; Elastase specific inhibitor; ESI; Peptidase inhibitor 3; PI 3; PI3; Protease inhibitor WAP3; SKALP; WAP3; WFDC14

To place an order, please [Click HERE](#).

Catalog #: A5441

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5266

Isotype: IgG

Swiss Prot: P19957

Purity: Affinity purification

For research use only.