

TPT1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP RIP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200 IP 1:20 - 1:50

RIP 1:20 - 1:50

Calculated MW: 20kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TPT1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bj

Synonym:

HRF; p02; p23; TCTP;

To place an order, please [Click HERE](#).

Catalog #: A5442

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7178

Isotype: IgG

Swiss Prot: P13693

Purity: Affinity purification

For research use only.