

CBR1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 30kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CBR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CBR; hCBR1; SDR21C1

Catalog #: A5446

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 873

Isotype: IgG

Swiss Prot: P16152

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the short-chain dehydrogenases/reductases (SDR) family, which function as NADPH-dependent oxidoreductases having wide specificity for carbonyl compounds, such as quinones, prostaglandins, and various xenobiotics. Alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).