

CPE

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 53kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CPE

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Carboxypeptidase E; Carboxypeptidase H; CPE; CPH; Enkephalin convertase;

Catalog #: A5458

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1363

Isotype: IgG

Swiss Prot: P16870

Purity: Affinity purification

For research use only.

Background:

This gene encodes a carboxypeptidase that cleaves C-terminal amino acid residues and is involved in the biosynthesis of peptide hormones and neurotransmitters, including insulin. It is a peripheral membrane protein. The protein specifically binds regulated secretory pathway proteins, including prohormones, but not constitutively secreted proteins. Mutations in this gene are implicated in type II diabetes.

To place an order, please [Click HERE](#).