

GGH

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 36kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GGH

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GH;

Catalog #: A5464

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8836

Isotype: IgG

Swiss Prot: Q92820

Purity: Affinity purification

For research use only.

Background:

This gene catalyzes the hydrolysis of folypoly-gamma-glutamates and antifolypoly-gamma-glutamates by the removal of gamma-linked polyglutamates and glutamate.

To place an order, please [Click HERE](#).