

HIF1AN

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP ChIP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200 IP 1:50 - 1:200

ChIP 1:20 - 1:100

Calculated MW: 40kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human HIF1AN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

FIH1;

To place an order, please [Click HERE](#).

Catalog #: A5466

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55662

Isotype: IgG

Swiss Prot: Q9NWT6

Purity: Affinity purification

For research use only.