

PPP2R5C

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 61kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PPP2R5C

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

B56G; PR61G;

Catalog #: A5480

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5527

Isotype: IgG

Swiss Prot: Q13362

Purity: Affinity purification

For research use only.

Background:

The product of this gene belongs to the phosphatase 2A regulatory subunit B family. Protein phosphatase 2A is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. The B regulatory subunit might modulate substrate selectivity and catalytic activity. This gene encodes a gamma isoform of the regulatory subunit B56 subfamily. Alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).