

VIL1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:93kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human VIL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

VIL; D2S1471;

Catalog #:A5494

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7429

Isotype:IgG

Swiss Prot:P09327

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of a family of calcium-regulated actin-binding proteins. This protein represents a dominant part of the brush border cytoskeleton which functions in the capping, severing, and bundling of actin filaments. Two mRNAs of 2.7 kb and 3.5 kb have been observed; they result from utilization of alternate poly-adenylation signals present in the terminal exon.

To place an order, please [Click HERE](#).