

DNAJB1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 38kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human DNAJB1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Hdj1; Sis1; HSPF1; Hsp40; RSPH16B;

To place an order, please [Click HERE](#).

Catalog #: A5504

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3337

Isotype: IgG

Swiss Prot: P25685

Purity: Affinity purification

For research use only.