

CDC25A

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 70kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human CDC25A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CDC25M3; MPIP1;

Catalog #: A5527

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 993

Isotype: IgG

Swiss Prot: P30304

Purity: Affinity purification

For research use only.

Background:

Tyrosine protein phosphatase which functions as a dosage-dependent inducer of mitotic progression. Directly dephosphorylates CDC2 and stimulates its kinase activity. Also dephosphorylates CDK2 in complex with cyclin E, in vitro.

To place an order, please [Click HERE](#).