

JAK2

Reactivity: Human Rat

Tested applications: IHC

Recommended Dilution: IHC 1:50 - 1:100

Calculated MW: 131kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human JAK2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bo

Synonym:

JTK10; THCYT3;

Catalog #: A5537

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3717

Isotype: IgG

Swiss Prot: O60674

Purity: Affinity purification

For research use only.

Background:

This gene product is a protein tyrosine kinase involved in a specific subset of cytokine receptor signaling pathways. It has been found to be constitutively associated with the prolactin receptor and is required for responses to gamma interferon. Mice that do not express an active protein for this gene exhibit embryonic lethality associated with the absence of definitive erythropoiesis.

To place an order, please [Click HERE](#).