

TFCP2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 57kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TFCP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

LSF; SEF; LBP1C; LSF1D; TFCP2C;

Catalog #: A5555

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7024

Isotype: IgG

Swiss Prot: Q12800

Purity: Affinity purification

For research use only.

Background:

This gene encodes a transcription factor that binds the alpha-globin promoter and activates transcription of the alpha-globin gene. The encoded protein regulates erythroid gene expression, plays a role in the transcriptional switch of globin gene promoters, and it activates many other cellular and viral gene promoters. The gene product interacts with certain inflammatory response factors, and polymorphisms of this gene may be involved in the pathogenesis of Alzheimer's disease.

To place an order, please [Click HERE](#).