

SH2B1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:50

Calculated MW: 79kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SH2B1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PSM; SH2B;

Catalog #: A5564

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 25970

Isotype: IgG

Swiss Prot: Q9NRF2

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the SH2-domain containing mediators family. The encoded protein mediates activation of various kinases and may function in cytokine and growth factor receptor signaling and cellular transformation. Alternatively spliced transcript variants have been described.

To place an order, please [Click HERE](#).