

RABL2A

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 26kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant Protein of human RABL2A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Rabl2a; AW549303; 1110031N17Rik;

To place an order, please [Click HERE](#).

Catalog #: A5575

Antibody Type:

Polyclonal Antibody

Species: Guinea pig

Gene ID: 11159

Isotype: IgG

Swiss Prot: Q9UBK7

Purity: Affinity purification

For research use only.