

WHSC1L1

Reactivity:Human Mouse

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:161kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human WHSC1L1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

NSD3; pp14328;

Catalog #:A5577

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:54904

Isotype:IgG

Swiss Prot:Q9BZ95

Purity:Affinity purification

For research use only.

Background:

This gene is related to the Wolf-Hirschhorn syndrome candidate-1 gene and encodes a protein with PWWP (proline-tryptophan-tryptophan-proline) domains. The function of the protein has not been determined. Two alternatively spliced variants have been described.

To place an order, please [Click HERE](#).