

## MLKL

---

**Reactivity:** Human

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:100

**Calculated MW:** 54kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human MLKL

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

9130019I15Rik;

*To place an order, please [Click HERE](#).*

**Catalog #:** A5579

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 197259

**Isotype:** IgG

**Swiss Prot:** Q9D2Y4

**Purity:** Affinity purification

For research use only.