

DHX9

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:1000 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 141kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human DHX9

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

LKP; RHA; DDX9; NDH2; NDHII;

Catalog #: A5585

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1660

Isotype: IgG

Swiss Prot: Q08211

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the DEAH-containing family of RNA helicases. The encoded protein is an enzyme that catalyzes the ATP-dependent unwinding of double-stranded RNA and DNA-RNA complexes. This protein localizes to both the nucleus and the cytoplasm and functions as a transcriptional regulator. This protein may also be involved in the expression and nuclear export of retroviral RNAs. Alternate splicing results in multiple transcript variants. Pseudogenes of this gene are found on chromosomes 11 and 13.

To place an order, please [Click HERE](#).