

DCLRE1C

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 78kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human DCLRE1C

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

e

Synonym:

SCIDA; SNM1C; A-SCID; hSNM1C; RS-SCID; DCLREC1C;

Catalog #: A5615

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64421

Isotype: IgG

Swiss Prot: Q96SD1

Purity: Affinity purification

For research use only.

Background:

This gene encodes a nuclear protein that is involved in V(D)J recombination and DNA repair. The protein has single-strand-specific 5'-3' exonuclease activity; it also exhibits endonuclease activity on 5' and 3' overhangs and hairpins when complexed with protein kinase, DNA-activated, catalytic polypeptide. Mutations in this gene cause Athabaskan-type severe combined immunodeficiency (SCIDA).

To place an order, please [Click HERE](#).